



LAMPIRAN KELENGKAPAN A1 SICS SET

KEMENKES RI AKL 11203220070

No.	Deskripsi	Tipe / Kode
1	Barraquer Eye Speculum, wire, 14mm open blades, length 4.0cm	ES-0122
2	Lieberman Eye Speculum, adj. mechanism, 15.0mm open wire blades	ES-0340
3	Graefe Muscle Hook, small, size 1	ES-0790
4	Sinskey II Lens Manipulating Hook,	HS-0690
5	Lewis Lens Loop, small	HS-1300
6	Barraquer Iris Spatula, round, smooth,	HS-0470
7	Nagahara Phaco Chopper, 1.25mm tip, offset blade, angled 60	HS-0360
8	Backhaus Towel Clamp, length 9.0cm	F-4690
9	Hartmann Mosquito Forceps, curved, length 9.5cm	F-4650
10	Hartmann Mosquito Forceps, straight, length 9.5cm	F-4640
11	Castroviejo Caliper, 20.0mm straight	M-0010
12	Iris Forceps, straight, with 1x2 teeth, 0.7mm, 90	F-4000
13	Iris Forceps, serrated, curved, 0.8mm, length 10.0cm	F-3930
14	McPherson Corneal Forceps, straight, with 5.0mm tying platform, 1x2 teeth, 0.12mm, 90	F-1040
15	McPherson Tying Forcepsc, straight, with 5.0mm tying platform, length 8.5cm	F-0950
16	McPherson Tying Forcepsc, angled, with 5.0mm tying platform, length 8.5cm	F-0960
17	Utrata Capsulorhexis Forceps, angled, triangular shaped tips, length 10.5cm	F-1820
18	Remky Superior Rectus Forceps, length 11.0cm	F-3780
19	Corneal Colibri Forceps, with tying platform 1x2 teeth, 0.12mm, 90	F-0230
20	Castroviejo Needle Holder, curved, 10mm x 1.0mm, smooth jaw, without lock, length 14.0mm	N-0120

21	McPherson Needle Holder, curved, 10.0mm x 0.5mm, smooth jaw, without lock, length 10.0cm	N-0080
22	Rainin Injection Cannula, angled, 7.0mm, 27G	L-0740
23	Hydrodissection Cannula, angled, 8.0mm, flat, front opening 25G	L-0750
24	Simcoe Irrigating/ Aspirating Cannula, 23G	L-1410
25	Knolle-Pearce Irrigating Vectis Cannula, with three front opening ports, 20G Loop	L-1090
26	Vannas Capsulotomy Scissors, curved, sharp tips, 5.0mm, length 8.5cm	S-0440
27	Castroviejo Universal Corneal Scissors, curved, blunt tips & small blades, cutting length 6.0mm, length 10.0cm	S-0160
28	Westcott Tenotomy Scissors, curved, blunt tips & medium blades, cutting length 8.00mm, length 11.5cm	S-0210