

# TRYPAN BLUE SOLUTION

5 x 1 ML Vial

## D-Blue

STERILE SOLUTION



**WARNING :** To be retail on the prescription of a Registered Medical Practitioner only

**NOTE :** For Single use only. Not for Injection. For intraocular use only  
Keep out of reach of children. Store in a cool dark place.  
Preservative free.

aadPri



### Composition :

Trypan Blue.....0.06% w/v  
Buffered Aqueous Vehicle.....q.s.



8 904278 700056

Packing: Packed in 1 ml Vial

SR.NO.	TEST	SPECIFICATION
1	Description	Dark blue colour solution filled in glass vial
2	Identification	Give positive test under assay
3	pH	Between 6.8 and 7.8
4	Extractable Volume	Between 1.0 ml to 1.15 ml
5	Sterility	Sterile
6	Assay	Between 0.054% w/v to 0.066% w/v of Trypan blue

AADPRI OPHTHALMIC COMPANY  
G-15, Mahavir krupa complex, Near C.P. Nagar Cross road,  
Bhuyangdev char.rasta, Ahmedabad- 380061, Gujarat, India.  
Mobile no: +91 9374862303  
Email: aadpri@gmail.com Web: www.aadpri.com

aadPri



EVERY EYE NEED EYE