





a1 MEDICAL PHACO SET

No.	CODE	DESCRIPTION	PICTURE
1	ES-0122	BARRAQUER Eye Speculum Wire, 14mm Open Blades, Length 4.0cm	
2	ES-0350	LIEBERMANN Eye Speculum, Adjustable Mechanism, 15.0mm V-Shaped Wire Blades	
3	HS.0690	SINSKEY II Lens Manipulator Hook, Ø 0.2mm Blunt Tip, Angled	
4	HS-0660	Lester IOL-Manipulator, Ø 0.25mm, 10.0mm Angled	
5	HS-0630	BECHERT Nucleus Rotator, Blunt Forked Tip, Y-Horizontal 7.0mm, Angled	
6	HS-1100	CHANG Phaco Chopper, a Combination of a Modified MicroFinger and a Pointed, Quick Chop Tip	
7	HS-0360	NAGAHARA Phaco Chopper, 1.25mm Tip, Offset Blade, Angled 60°	
8	HS-1120	Culler Iris Spatula 1.0mm Wide	
9	F-1040	McPHERSON Corneal Forceps, Straight, with 5.0mm Tying Platform 1x2 Teeth, 0.12mm, 90° Angle, Length 8.5cm	



10	F-2470	McPHERSON Suture Forceps, Straight, with 5.0mm Tying Platform, Length 10.0cm	
11	F-1820	UTRATA Capsulorhexis Forceps, Angled, Trianglurar Shaped Tips, Length 10.5cm	
12	F-0980	Bechert-McPHERSON Tying Forceps 10.0mm, Angled, with 5.00mm Tying Platform, Length 8.5cm	
13	F-0230	Corneal Colibri Forceps with Tying Platform, 1x2 Teeth 0.12mm, 90° Angled, Length 7.5cm	
14	N-0120	CASTROVIEJO Needle Holder Curved 14.0mm	
15	N-0080	McPherson Needle Holder Curved 10.0mm X 0.5mm, Smooth Jaw, Without Lock, 10.0cm	
16	S-0490	Gills-Vannas Capsulotomy Scissors, Angled Forward, Sharp Tips, 7.0mm Blade, Length 8.5cm	
17	S-0010	Castroviejo Universal Corneal Scissors, curved, blunt tips & small blades, cutting length 5.0mm, length 11.5cm	
18	L-0740	Rainin Injection Cannula, angled, 7.0mm, 27G	
19	L-0750	Hydrodissection Cannula, angled, 8.0mm, flat, front opening, 25G	
20	L-1410	SIMCOE Irrigating / Aspirating Cannula 23g	

