






SL-1S Optics(3X)



SL-2S Optics(5X)



Clear and Stereoscopic View

❖ Parallel type, Galilean system

Superior optic quality, provides stereoscopic and comfortable view. With broad beam, wide field and high depth of field

❖ Warm LED Illumination, Better Conformity, Longer Life Time

Work 50 thousand hours, and suitable for human tissue.

❖ Precision Machining Makes Smooth and Flexible Movement

❖ Brightness Control Continuously

❖ Interface for laser and applanation tonometer reserved

❖ Exquisite shape design, compact mobile base



See More, See Future!



Technical Parameter

Items	Model	SL-1S	SL-2S
Optical Part			
Microscope Type	Galilean Type System		
Magnification Changer	Revolving Drum		
Magnification	10X,16X,25X	6X,10X,16X,25X,40X	
Eyepiece Magnification	12.5x		
Visual Angle of Eyepiece	13°		
PD Range	52mm-85mm		
Diopter	±7D		
Visual Field Diameter	1 0 X (2 2 m m) , 1 6 X (1 4 m m) , 25X(8.5mm)	6 X (3 3 m m) , 1 0 X (2 2 m m) , 1 6 X (1 4 m m) , 25X(8.5mm), 40X(5.5mm)	
Illumination Part			
Illumination Type	Tower Type		
Light Source	Warm LED		
Brightless control	Continuously adjustable		
Illumination tilting	5°,10°,15°,20°		
Slit Width	0mm-14mm continuously adjustable		
Slit Height	1mm-14mm continuously adjustable		
Slit Angle	0°~180°continuously adjustable		
Slit Aperture Size	Φ0.2,Φ1,Φ3,Φ5,Φ10,Φ14,mm)		
Filter	Heat Aborption, Grey, Redfree(Green), Cobalt Blue		
* Voltage: 100~240V 48~62Hz Compatible			

Main MediTech Co., Ltd

📍 Add: 9#&10# of No. 776, JingWei Avenue, YuZhong District, ChongQing, China

🌐 Website: www.mainmeditech.com

✉ Email: info@mainmeditech.com

☎ Tel: +86 23 6373 2637

☎ Supervisor Tel : +86 13618270311

See More, See Future!